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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION

The massage pillow is used to massage the neck, spine, back and legs. It provides proper muscle relaxation, giving a feeling of comfort and relaxation. It helps with back pain, chronic muscle tension, varicose veins, and also relaxes the muscles of the neck, shoulders or back. It warms up the muscles, improves circulation, as well as overall physical and mental condition. The additional function of heat distribution breaks down compacted muscles due to stress, incorrect posture or tension. It has a positive effect on the nervous system, calms, relaxes and unwinds. Its small size means that you can take it with you on a trip, to work, and even use it in the car thanks to the included car cable, the pillow has a strap on the back that allows it to be mounted to the headrest in the car.

INSTRUCTIONS FOR USE:

- Remove packaging from product.
- Place the pillow in the desired location: behind the neck, behind the back, under the lower back or under the legs.
- Sit or lie down comfortably so that the massage heads adhere well to the part of the body being massaged.
- Plug the massager into a wall outlet.
- Press the power button.
- The device will automatically activate the heated massage function.
- You can change the rotation direction of the massage heads by pressing the function button.
- The massage usually lasts about 15 minutes.
- Once you have finished massaging, you can turn off the device manually.
- Unplug the power cord from the wall outlet or cigarette lighter socket.
- Store the pillow in a dry and safe place, protecting it from dust and moisture.

TIP

- Using the device for purposes other than its intended purpose and/or using the device in any other way may involve various risks.
- Use the device only for its intended purposes.
- The procedures described in this operating manual must be followed. Claims of any kind resulting from damage caused by improper use, unauthorised repairs, unauthorised changes or the use of unauthorised spare parts are excluded. The risk lies solely with the user.
- If you experience any skin discomfort or irritation, discontinue use immediately.
- Do not pull or excessively stretch the cord while using the appliance.
- The device is not intended for children - do not allow them to play with the massager.

TECHNICAL DATA

- 4 rotating massage heads
- Removable pillowcase
- Car charger: 2A
- Power: 24 W
- Voltage: AC 110-240V
- Two modes: normal and infrared
- Dimensions: 30cm x 19cm x 10cm
- Car charger
- Mains charger

SAFETY INSTRUCTIONS

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Never immerse the device in water.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.
- The product is not intended for children to play with. Children must not play with the device.
- Cleaning and user maintenance shall not be performed by children without supervision.
- The product itself and its accessories are not toys, please keep them out of the reach of children to avoid danger.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.

- It is forbidden to use the product if any of its parts are damaged. In the event of damage to the cable, it is forbidden to make repairs yourself.
- Do not use if the power cord or power adapter is damaged or if it is not properly plugged into the outlet.
- Do not disassemble the device yourself.
- Do not hit or drop this product. Stop using it if it is damaged in a fall or other mechanical damage.
- Store the device in a dry place.
- Before using the device, make sure that the mains voltage matches the required operating voltage of the device.
- Always use the product as intended.

ELECTRIC SHOCK HAZARD

- Connect the device only to a properly installed and grounded power outlet. The mains voltage must match the data on the nameplate of the device / mains adapter.
- Make sure that the power cable does not get wet or damp during operation. The cable must be routed so that it is not pinched or damaged.
- Keep the power cable away from hot surfaces.
- Do not perform any repairs to the device yourself.
- Any repairs may only be performed by service or by qualified specialists.
- Before cleaning the device, unplug the power supply from the wall outlet.
- Do not use the device in a damp or wet environment, such as a bathroom.
- Have a damaged plug or power cable repaired immediately by a qualified specialist or service center in order to avoid any danger.
- Do not insert any needles or sharp objects into the device.

WARNING! INJURY RISK!

- This device must not be used to massage infants, small children, helpless persons or animals.
- Never use this device if you experience dislocation or pain in the area being massaged.
- Stop the massage immediately if it becomes uncomfortable or causes pain.
- Do not use the device on swollen, burned, irritated, diseased or injured areas of the skin or body.
- If you have any doubts regarding use, please consult your doctor.
- Do not fall asleep while using your device.
- Pregnant women, people with pacemakers, heart disease, diabetes, neurological or skin problems should consult their physician before using the device.
- Do not fall asleep during the massage – the device works mechanically and requires control.

WARNING – MATERIAL DAMAGE!

Use the appliance indoors only.

Only use the device with the supplied mains adapter.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the appropriate office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Electric shock warning! Danger to life!

We reserve the right to make changes to the text, design and technical data of the product without notice.

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